# ACE AFM-5

Alkoholtester Breathalyzer

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### Author's Liability:

All information herein was created conscientiously and under the assumption of correctness. Hence, the author cannot be held liable for the use of any kind whatsoever and/or any ensuing consequences. Furthermore, the author is not liable for any costs or expenses for damage that result from any information in this document. Any information is subject to change without prior notice.

### Introduction:

Thank you for purchasing our ACE AFM-5 Breathalyzer. It measures alcohol concentration in human breath using an electrochemical sensor and calculates the blood alcohol level. The test results are very reliable and have a tolerance of  $\pm 0.050\%$  (per mill) with reference to  $1.0\%$ . The device must only be used pursuant to the instructions mentioned in this manual.

### Short description:


The ACE AFM-5 Breathalyzer is an easy-to-use alcohol tester with an electrochemical sensor. It is used in combination with a smartphone (only iOS and Android system possible / no tablets). Before performing the test you have to download the "Alcofind" app from the App Store (iOS) or from Google Play (Android) and install the app on your smartphone. The test result can be interpreted 1:1 (so if the shown test result is, e.g., 0.360, then this means 0.360 per mill).  $\%$  (per mill) = measurement unit in the EU, in the App "Alcofind" you can change the measurement unit into % BAC (USA) or into mg/L (UK).

The fuel cell sensor inside the ACE AFM-5 Breathalyzer can be calibrated. For a calibration service you can send the device to ACE Instruments (see [www.ace-instruments.net](http://www.ace-instruments.net)). For each breath test a new mouthpiece is required.

### Scope of delivery:

- 1 Digital fuel cell, ACE Instruments AFM-5 Breathalyzer
- 1 AAA battery (1.5 V)
- 6 Mouthpieces
- 1 Carrying case
- 1 Quick Start Guide (German & English)
- 1 Manual (German & English)

### Label of individual components:

This way you receive results that are as accurate as possible:

- a) Wait at least 15 minutes after your last sip of alcohol. Do not eat, drink or smoke during that time.
- b) Do not open the casing.
- c) Do not keep the ACE AFM-5 Breathalyzer in cold, hot, moist or badly ventilated places.
- d) Keep the ACE AFM-5 Breathalyzer out of the reach of children.
- e) Conduct three inlet tests after long storage before conducting the actual alcohol test.
- f) Do not drop the ACE AFM-5 Breathalyzer.
- g) Change the battery if it is depleted. Diminished battery capacity leads to incorrect test results.
- h) Do not conduct tests in loud (sound waves disrupt the sensor), drafty or poorly ventilated places.
- i) After a test result of 1.000 per mill (per mill = EU setting) / 0.100 % BAC (% BAC = US setting) / 0.500 mg/L (mg/L = UK setting) or higher the sensor needs rest time. In this case wait at least three minutes before starting the next test.
- j) Do not conduct tests in very cold, very hot or an unventilated, smoky ambience.
- k) Store the ACE AFM-5 Breathalyzer in a way to avoid penetration by dust or liquids.
- l) Blowing too hard may cause damage to the sensor as saliva flows into it.

## Preparation for the test:

### 1. Installing a battery

**Step 1.** Push up the bottom part of battery cover to open it.

**Step 2.** Install an AAA battery in the battery compartment while taking note of the battery polarity.

**Step 3.** Insert bulging parts of the battery cover into the slits, and then push down its cover.

### 2. Inserting a mouthpiece

Insert a mouthpiece into the slot on the side of the alcohol tester. Ensure whether it securely fits or not. It is highly recommended to use an extra mouthpiece for every test for sanitary reasons.

### 3. Installing the ACE Instruments AFM-5 app

**Step 1.** Download the “**Alcofind**” app from App store (iOS) or Google Play Store (Android). Note the minimum smartphone requirements: Either an iPhone 4S running iOS 7 and above or Android devices running OS 4.3 and above.

**Step 2.** Open the app. Select language and unit, then press “**OK**” button. **NOTE:** For changing language and unit later, go to setting at the bottom menu.

### Operation with App (App mode):

**Step 1.** Open the app. **NOTE:** Bluetooth must first be switched on in order to connect to the ACE Instruments AFM-5. (Go to the settings on your smartphone and turn on Bluetooth.)

**Step 2.** To turn on the ACE Instruments AFM-5, press and hold Power button for a second until the Blue, Green and Red LEDs appear on the device as shown in the picture below.

**Step 3.** The Blue LED will blink during pairing process. If pairing is successful, Green and Red LEDs will appear on the device and the screen on your smartphone will show “**START TEST**”. Tap “**START TEST**” on your smartphone (or press Power button on the ACE Instruments AFM-5).

**NOTE:** In case of the first pairing process, a PIN Code request will be displayed. Type in “**000000**” as PIN Code.

**NOTE:** The app is set to remember and search the device which has paired last time for 30 seconds. If the app is not able to find the previously paired ACE Instruments AFM-5 for 30 seconds, the app will show a request message. If you want the app to find another ACE Instruments AFM-5, press “**OK**”, and if you want the app to keep finding the previously paired ACE Instruments AFM-5, press “**Cancel**”.

**Step 4.** Choose a user. Each test result will be saved separately according to the selection of the user. All history of saved test results are viewable from the Graphic Icon at the bottom of the menu.

**Step 5.** The circular bar on the app indicates the progress of

warm-up. When warming up gage reaches 95, take a deep breath and get ready for blowing.

**Step 6.** When “**Start blowing**” is displayed, blow for 5 seconds into the mouthpiece until the gage reaches 100 or you feel the air has been pulled through the air pump in the ACE Instruments AFM-5 with a click.

**NOTE:** If the camera feature is enabled, the camera will automatically capture the image of a user once the blowing has been completed.

**Step 7.** If a sufficient breath sample is taken correctly, the “**Analyzing data**” will be shown on the app. **NOTE:** If insufficient sample is taken, an error message will be displayed. To test again, press “**OK**” on the app. **NOTE:** In the graph, the test results performed without app will be shown as “**D**”.

**Step 8.** After a few seconds the estimated blood alcohol concentration will be displayed on the app in per mill (3 digits, 2 decimals/per mill = default setting). The result may vary depending on the concentration which has been measured during the last test. In this page, the app enables you to make a memo which is viewable with test results from the Graph icon and the test results can be also shareable via Text message, Facebook, or Twitter.

**Step 9.** Tap either “**Retest**” or “**BAC level**” on the app (or press Power button on the ACE Instruments AFM-5) to test again without selecting a user. If you want to start test again from the beginning, tap the Home Icon at the bottom menu.

**Step 10.** Press and hold Power button on the ACE Instruments AFM-5 for 2 seconds to turn it off when you have finished testing.

**NOTE:** The device will turn off automatically if the device is not in use for 2 minutes. If you exit the app while the device is working with the app, all Blue, Red and Green LEDs on the device will flicker together. Press Power Button on the device to restart.

### Operation without App (Standalone mode):

**NOTE:** To activate Standalone mode, you have to pair the ACE Instruments AFM-5 with the app and set “**Limit Line**” first. (Pair with the app and go to the setting on the app to set “**Limit Line**”.)

Once “**Limit Line**” is set, the set value will be sent to the ACE Instruments AFM-5, and the test result on Standalone mode will be displayed with indications of Red or Green LEDs according to the set value of “**Limit Line**”. (After setting the value of “**Limit Line**”, you should move into Home by tapping the Home Icon at the bottom menu on the app in order to save and transfer the set value to the ACE Instruments AFM-5.)

**NOTE:** If the battery is low or the device is in unacceptable temperature range, all Blue, Red and Green LEDs will flicker together. Replace a new battery or try again after the device has reached a temperature within the operating range (5-40 °C).

**Step 1.** To turn on the ACE Instruments AFM-5, press and hold Power button for a second until the Blue, Green and Red LEDs appear on the device.

**Step 2.** If the app on your smartphone is off and there is no pairing signal from the app, the device will move into the Standalone Mode and starts the warm-up automatically. A Blue LED keeps blinking until the warm-up process is completed. **NOTE:** If the device doesn't have the set value of “**Limit line**”, all Blue, Red and Green LEDs will flicker together after the warm-up process. Pair the device with the app and set “**Limit Line**” on the setting menu of the app. If you turn on the app on your smartphone while using the device in the Standalone Mode, the device will move into the App Mode and start pairing automatically. (If you turn on the app while the device is under the blowing or analyzing process in the Standalone Mode, the device will move into the App Mode after completing these processes.)


**Step 3.** When the Blue LED stops blinking and remains on, it means the device is ready for blowing. Take a deep breath and blow for 5 seconds into the mouthpiece until the Blue LED disappears.

**Step 4.** If a sufficient breath sample is taken correctly, Green and Red LEDs will blink in turns and the test result will be displayed in a few seconds. **NOTE:** If an insufficient breath sample is taken, blue, red and green LEDs will flicker together. Press Power button to retest.

**Step 5.** Test results will be displayed with indications of Red or Green LEDs according to the set value of “**Limit Line**” as follows. Press Power Button to retest.

## Statement of measurement levels:

Example: If you have set 0.5 per mill in the limit line of the App and your current breath sample is 0.3 per mill then the green light blinks at an interval of 1 second.

\*e.g. 0.5‰, per mill = default measurement, can be changed to %BAC OR mg/l

**NOTE:** If "**Limit Line**" is set as 0.20‰, the indication of pass will be displayed as 2 levels - PASS (HIGH) & PASS (MIDDLE).  
If "**Limit Line**" is set as 0.10‰, the indication of pass will be displayed as 1 level - PASS (HIGH).

**NOTE:** The ACE Instruments AFM-5 saves up to 30 test results and the oldest test results will be deleted automatically if memory is full. Once the ACE Instruments AFM-5 pairs with the app, the saved test results will be sent to the app, and at the same time the saved test results will be automatically deleted from the ACE Instruments AFM-5. All transmitted test results are viewable from the Data Icon at the bottom of the app menu. (In the graph, the transmitted test results will be shown as "D").

**Step 5.** Press Power button on the ACE Instruments AFM-5 to retest.

**Step 6.** Press and hold Power button for 2 seconds to turn it off when you have finished testing. **NOTE:** The device will turn off automatically if the device is not in use for 2 minutes.

## ACE Instruments recommends:

- Always use a new mouthpiece for each person and test that is

conducted. ACE Instruments assumes no liability for any damages which arise when several people use the same mouthpiece for multiple tests.

- Only use the provided single packed and hygienic mouthpieces - more are available at [www.ace-instruments.net](http://www.ace-instruments.net). During each test the ACE AFM-5 Breathalyzer's internal and external parts come into contact with human secretion (saliva, blood, sweat, etc.)!
- Use disposable gloves when pulling out a mouthpiece used by another person! You can buy appropriate gloves at ACE Technik ([www.ace-technik.de](http://www.ace-technik.de)).

## Important notes:

**When you have eaten, drunken alcohol or smoked a cigarette, wait at least 20 minutes before starting the measurement.** If the test result is below 0.050‰ (per mill) OR 0.005% BAC OR 0.025 mg/l (depending on the setting) then the device shows 0.000‰ OR 0.000% BAC OR 0.000 mg/l on the display. This happens because at a low range of alcohol concentration, and thus possible to measure other materials in your mouth.

## Only the usage of original ACE mouthpieces...

- Guarantees test result accuracy
- Guarantees 100 % hygieny
- Guarantees full warranty of this device

## ACE Instruments advises:

Even if you get a result meaning "**safe range**", you should not drive after consuming alcohol or do any other activities that are severely impaired by alcohol consumption.

## Notes to the battery law (BattG):


The device associated with this operating manual contains one AAA Alkaline battery (1.5 V) in delivery condition. This battery is not made for recharging. The use of a rechargeable battery is not possible with this device and can damage it. If the included battery is depleted, it must not be thrown into the household

waste. Used batteries may contain harmful substances that can cause environmental or health damage. Please return all depleted batteries to electrical appliance stores, who collect them for free.

### Battery exchange:

To replace the battery, push the battery cover on the rear of the unit up and replace the discharged battery with a new one. Be aware of the polarity when inserting the battery (+/- pole). Put the battery compartment cover on the rear panel and press it gently on the device until it clicks. In the EU the return of used batteries is mandatory by law. All depleted batteries and also rechargeable batteries will be recycled. This way valuable resources like iron, zinc and nickel can be retrieved and are reusable then. Battery recycling is the easiest lesson in environmental protection – so please support it!

The shown garbage bins mean that batteries or rechargeable batteries must not be thrown into the household waste. The signs below the garbage bins are acronyms for the chemical substances inside the battery:

### Calibration:

To maintain the accuracy of the ACE AFM-5 Breathalyzer, you can send the device for a calibration service to ACE Instruments (see [www.ace-instruments.net](http://www.ace-instruments.net)). The calibration is not mandatory, but the accuracy of ACE AFM-5 Breathalyzer will decrease steadily without it. We recommend sending it as a parcel, so that in cases of loss an inquiry or proof of delivery is possible. Due to changes in the product, the procedures stated can change. This can take place without prior notice or adaptation of the operating manual and is no cause for complaint, return or price reduction. The user is responsible for recognizing discrepancies between the operating manual and the actual handling of the ACE Instruments AFM-5 Breathalyzer and must react accordingly.

### Technical data\*:

Sensor:	Fuel Cell sensor
Accuracy:	+/- 0.05‰ (per mill) with reference to 1.0‰ at 25°C (per mill = measurement unit in the EU) OR +/- 0.005% BAC (with reference to 0.1% BAC at 25°C (% BAC = measurement unit in the USA OR +/- 0.025mg/L with reference to 0.5mg/L at 25°C (mg/L = measurement unit in UK)
Detection Range:	0,000 - 4,000‰ OR 0,000 - 0,400% BAC OR 0.000 - 2.000 mg/L (per mill = measurement unit in the EU), in the App "Alcofind" you can change the measurement unit into % BAC (USA) or mg/L (UK)
Warm-up Time:	Within 15 seconds at 0.5‰, it may vary depending on the value of per mill / BAC which was measured
Connectivity:	Bluetooth 4.0 LE
Dimension:	68mm x 41.5mm x 16.9mm (L x W x D)
Weight:	43g (battery included)
Power Supply:	1 piece of 1.5V AAA size Alkaline battery
Working Temperature:	5°C - 40°C
Storage Temperature:	0°C - 40°C
Calibration:	Every 12 months or after 500 tests

\*The technical data are subject to change without prior notice for functional improvements.

## Warranty and liability:

### WARRANTY

ACE Instruments (ACE) provides a warranty of two (2) years after purchase for material and processing errors for the product (see attached warranty and liability declaration) in all versions and designations; and a warranty of one (1) year for the sensor. This warranty includes the electronics as well as the casing or sensor without sampling components. Components that are inevitably destroyed upon effectiveness (such as fuses) are not part of this warranty. The warranty is valid from the date of delivery or from the date of first activation (always the shorter period), if not stated otherwise in the accompanying ACE product literature.

### MISUSE

OTHER APPLICATIONS OF, OR MODIFICATIONS TO, THIS BREATHALYZER CONSTITUTE A MISUSE AND ENTAIL A RISK OF PERSONAL INJURY AND/OR DAMAGE TO THE DEVICE. ACE SHALL ASSUME NO LIABILITY FOR INJURY/DAMAGE CAUSED BY MISUSE OF THE DEVICE. THE DEVICE IS NOT INTENDED FOR COMMERCIAL USE.

ACE SHALL NOT BE HELD LIABLE FOR ANY DAMAGES THAT OCCUR AS A RESULT OF ANY UNAUTHORIZED MODIFICATIONS, REENGINEERINGS OR AMENDMENTS MADE TO THIS BREATHALYZER OR ITS IMPROPER USE.

## Limitation of warranty and liability:

### 1. PRODUCT DESCRIPTION

ACE REJECTS ANY FURTHER EXPRESSLY STATED OR IMPLIED WARRANTY, INCLUDING THE MARKETABILITY AND SUITABILITY FOR A DETERMINED PURPOSE. THIS WARRANTY ESPECIALLY EXCLUDES CALIBRATION. AFTER CHANGING THE SENSOR, CALIBRATION IS ABSOLUTELY NECESSARY.

IF AN EXCHANGE OF THE SENSOR OR A CALIBRATION IS REQUIRED FOR A PRODUCT, THE USER MAY RECEIVE A COST ESTIMATE. SHOULD THE USER REJECT THE ESTIMATE, ACE

WILL SEND THE PRODUCT BACK UNREPAIRED/UNCALIBRATED AND THE PRODUCT WILL BE LABELLED ACCORDINGLY.

SHOULD THE PRODUCT NOT FULFIL THE WARRANTY, THE CLAIM OF THE BUYER AND THE OBLIGATIONS OF ACE ARE LIMITED EXCLUSIVELY AND AT THE DISCRETION OF ACE TO AN EXCHANGE OR REPAIR OF THE AFFECTED GOODS OR TO A REIMBURSEMENT OF THE PURCHASING PRICE FOR THE FAULTY GOODS.

### 2. APPLICATION RESTRICTION

ACE IS ON NO ACCOUNT LIABLE FOR ADDITIONAL DAMAGES, INCIDENTAL DAMAGES OR CONSEQUENTIAL DAMAGES SUCH AS, AMONG OTHERS BUT NOT LIMITED TO, LOST PROFIT AND DOWNTIME THAT RESULT FROM SELLING, PRODUCING OR USING PRODUCTS, WHICH ARE PART OF THIS APPLICATION RESTRICTION. THIS IS INDEPENDENT OF ANY CLAIMS, SUCH AS FROM CONTRACTS, NEGLIGENCE OR OTHER CULPABLE ACTIONS.

THE PRODUCT, WHICH IS ACCOMPANIED BY THIS WARRANTY AND LIABILITY DECLARATION, NO MATTER IN WHICH VERSION OR DESIGNATION, CAN ESPECIALLY NOT BE USED LEGALLY BINDING TO DETERMINE CAPABILITY TO OPERATE VEHICLES OR TO EXECUTE ACTION THAT COULD LEAD TO DAMAGE TO ONESELF OR OTHERS UNDER THE INFLUENCE OF ALCOHOL.

THE THOROUGH CHECK FOR DAMAGE OF ALL PRODUCTS UPON RECEIPT, A CALIBRATION SUITABLE FOR THE GIVEN USE BY THE BUYER AND THE PROPER USE, MAINTENANCE AND SERVICE ACCORDING TO THE INSTRUCTIONS OF THE ACCOMPANYING PRODUCT LITERATURE BY ACE ARE MANDATORY PREREQUISITES FOR THIS WARRANTY. MAINTENANCE AND SERVICE WORK BY NON-QUALIFIED PERSONS, AS WELL AS USING NON-APPROVED CONSUMABLES OR SPARE PARTS ARE A MISUSE, AND ALL WARRANTY RIGHT SHALL SEIZE IMMEDIATELY. AS FOR ALL TECHNICALLY CHALLENGING PRODUCTS, IT IS A FUNDAMENTAL PREREQUISITE AND MANDATORY FOR THE WARRANTY BY ACE THAT ALL USERS OF THIS PRODUCT ARE FAMILIAR WITH THE PROPERTIES AND RESTRICTIONS OF THE PRODUCT AS STATED IN THE RESPECTIVE LITERATURE.

## DISCLAIMER

TO THE EXTENT PERMITTED BY LAW, THIS WARRANTY AND THE REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, REMEDIES AND CONDITIONS, WHETHER ORAL OR WRITTEN, STATUTORY, EXPRESS OR IMPLIED. AS PERMITTED BY APPLICABLE LAW, ACE SPECIFICALLY DISCLAIMS ANY AND ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND WARRANTIES AGAINST HIDDEN OR LATENT DEFECTS. IF ACE CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES BE LIMITED IN DURATION TO THE DURATION OF THE EXPRESS WARRANTY AND TO THE REPAIR OR REPLACEMENT SERVICE AS STATED ABOVE. EXCEPT AS PROVIDED IN THIS WARRANTY AND TO THE MAXIMUM EXTENT PERMITTED BY LAW, ACE IS NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ANY LEGAL THEORY.

Some countries, states and provinces do not allow the exclusion or limitation of incidental or consequential damages or allow limitations on how long an implied warranty or condition may last, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary by country, state or province. This Limited Warranty is governed by and construed under the laws of the country in which the product purchase took place.

### 3. ADVICE RESTRICTION

The buyer confirms that he has checked the intended use and suitability of purchased goods for any purpose in sole responsibility. Between the two parties it is expressly agreed upon that all technical or other advice by ACE regarding use of the goods or services is free of charge and was made with to the exclusive risk of the buyer. Therefore no obligation or liability can arise to ACE from any recommendation or their result.

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ACE Handels und Entwicklungs GmbH